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Sinuscan 201

Sinuscan 201

Inhaltsangabe

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12. EMC-Information

The Sinuscan 201 meets the requirements of the EMC-standard IEC 60601-1-2 for Medical electrical equipment. Medical electrical equipment needs special precautions regarding EMC and need to be installed and put into service according to the EMC information provided here.

Fixed RF transmitters, portable and mobile RF communications equipment can affect to Sinuscan 201 and the tables 3 and 4 are guiding to prevent from interferences.

The device is suitable for use in all establishments, including domestic establishment and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes

Table 1: Emission	21
Table 2: Immunity	21
Table 3: Immunity in RF-field.....	22
Table 4: Separation distances	22

Table 1: Emission

Guidance and manufacturer's declaration – electromagnetic emissions		
The SINUSCAN 201 is intended for use in the electromagnetic environment specified below. The customer or the user of the type device should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishment, including domestic establishment and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes
Harmonic emissions IEC 61000-3-2, Class D	Not applicable, active input power <50 W	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

Table 2: Immunity

Guidance and manufacturer's declaration – electromagnetic immunity			
The SINUSCAN 201 is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	< ±1 kV for input/output lines	Occasional false echoes may occur even with mains power quality of a typical commercial or hospital environment
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_i (>95% dip in U_i) for 0.5 cycle 40 % U_i (60% dip in U_i) for 5 cycle 70 % U_i 30 % dip in U_i) for 25 cycle <5 % U_i (>95% dip in U_i) for 5 sec	<5 % U_i (>95% dip in U_i) for 0.5 cycle 40 % U_i (60% dip in U_i) for 5 cycle 70 % U_i 30 % dip in U_i) for 25 cycle <5 % U_i (>95% dip in U_i) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment

Note: U_i is the a.c. mains voltage prior to application of the test level

Table 3: Immunity in RF-field

Guidance and manufacturer's declaration – electromagnetic immunity			
The SINUSCAN 201 is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 V	Portable and mobile RF communications equipment should be used no closer to any part of the SINUSCAN 201, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1,2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2,3 \sqrt{P}$ 800 MHz to 2,5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
<p>a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the SINUSCAN 201 is used exceeds the applicable RF compliance level above, the SINUSCAN 201 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the SINUSCAN 201.</p>			

Table 4: Separation distances

Recommended separation distances between portable and mobile RF communications equipment and Sinuscan 201			
The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Sinuscan 201 as recommended below, according to the maximum output power of the communications equipment			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter		
	150 kHz to 80 MHz	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

13. Anhang

13.1 Bescheinigungen

QUALITY SYSTEM

CERTIFICATE

Directive 93/42/EEC

Manufacturer: Oriola Oy
Orionintie 5
02100 Espoo, Finland


Coverage of Certificate: Final inspection

Product category: Ultrasound diagnostic equipment


Valid until: 20th May 2007.

The manufacturer's quality system for the final inspection of the aforesaid product category has been evaluated and meets the provisions of Council Directive 93/42/EEC as set out in Annex VI Section 3. This approval is valid until the expiry date provided that the manufacturer fulfils the obligations imposed by Annex VI in Directive 93/42/EEC. This Certificate is based on decision no. TUO 1027-188-P0.

Tampere, 23rd June 2004

Markku Helminen

Kaarle Kylmä

Certificate no.
TUO 1027-188


Translated from Finnish original.

Notified Body no. 0537:
VTT Industrial Systems
P.O.Box 1306 (Tekniikankatu 1)
FIN-33101 TAMPERE
Tel.+358 3 316 3111

DECISION

on the changes to the approved quality system as specified in Annex VI, Section 3 of Council Directive 93/42/EEC concerning medical devices, or to the product range covered.

Decision no.:	TUO 1027-188-P1
Manufacturer:	Oriola Oy, Orionintie 5, FIN-02100 Espoo, Finland
Date of manufacturer's notification:	3 rd March 2005
Procedure and product category:	Final inspection of ultrasound diagnostic equipment.
Description of the changes:	Extension of the product category to cover the following product: - Sinuscan 201, class IIa
The decision:	The manufacturer's quality system for the the products above has been assessed and it meets the requirements in Annex VI of Medical Device Directive 93/42/EEC. The decision is based on the report no. NB-1027-MR02.
Certificate related to the decision:	AUT 1027-188
Validity:	This decision is valid until 20 th May 2007 unless the validity of the related certificate is changed.
Date:	Tampere, 11 th April 2005  Markku Helminen Kaarle Kylmä

VTT Industrial Systems is Notified Body no 0537 under Council Directive 93/42/EEC.

VTT INDUSTRIAL SYSTEMS

Tekniikankatu 1, Tampere
P.O.Box 1306, FIN-33101 TAMPERE
FINLAND

Tel. +358 3 316 3111
Fax +358 3 316 3365

name.surname@vtt.fi
www.vtt.fi
Business ID 0244679-4

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DECLARATION OF CONFORMITY

Manufacturers Name:

Oriola Oy
Orionintie 5
FIN-02101 ESPOO
FINLAND

Declares that the product:

Sinuscan ultrasound instrument, type 201, with power supplies
Friwo FW7333M/12 or Friwo FW7333M/12 UK or Friwo FW7555M/12
Conforms to the following European Union directives and standards
identified in this declaration.

EU Directives:

93/42/EEC, Medical devices


Standards:

ISO 14971:2003, Medical devices - Application of risk management to medical devices
IEC 60601-1+A1+A2, Medical electrical equipment - part 1: general requirements for safety
IEC 60601-1-1, General requirements for safety - Collateral standard: Safety requirements for medical electrical systems
IEC 60601-1-2, Medical electrical equipment: Part 1: General requirements for safety. 2. Collateral Standard: Electromagnetic compatibility - Requirements and tests
IEC 60601-2-37, Medical electrical equipment - Particular requirements for the safety of ultrasonic medical diagnostic and monitoring equipment
IEC 61157, Requirements for the declaration of the acoustic output of medical diagnostic ultrasonic equipment

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Mariachi Oy