



TS7000DV SURGICAL TABLE REDEFINING PATIENT POSITIONING IN **DA VINCI** SURGERY

The TS7000dV Surgical Table is:

- Uniquely compatible with the **da Vinci Xi** Surgical Robot, yet also serves as a versatile general surgical table for all specialties.
- Synchronized motion with the **da Vinci Xi** System for intraoperative re-positioning.
- Dramatic range of motion with high patient weight capacity in all table articulations.



INTEGRATED TABLE MOTION

Integrated Table Motion enables the **da Vinci Xi System** to connect to Baxter's TS7000dV Surgical Table so the patient can be dynamically positioned while the surgeon operates. This allows the surgical team to efficiently manage the patient and utilize gravity for optimal exposure and access to the target anatomy.

Iso centered motion: a combination of Trendelenburg or reverse Trendelenburg and slide.



DYNAMIC ACCESS & EXPOSURE

Use gravity to expose anatomy, maximize reach, and deliver tissue at an ideal working angle.



EFFICIENT WORKFLOW FOR EVERY CASE

Save time by eliminating the need to undock, resulting in shorter case times and the ability to provide intraoperative patient supine relief.



SEAMLESS INTEGRATION WITH PRECISE FEEDBACK

The **da Vinci Xi System** and TS7000dV Surgical Table connect wirelessly, and provide feedback on table angles for the surgical team.

For more information please contact your local Baxter representative.

hillrom.com

Markkinoija Suomessa



Mediq Suomi Oy

Vuoritontunkuja 6, 02200 Espoo

puh. 020 112 1510

asiakaspalvelu@mediq.com, www.mediq.fi

da Vinci and da Vinci Xi are trademarks or registered trademarks of Intuitive Surgical, Inc.

The TS7000dV Surgical Table is intended to be used by clinicians and medically qualified personnel. This medical device is a regulated healthcare product which, pursuant to such regulation bears a CE mark. Baxter recommends that you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices. Baxter reserves the right to make changes without notice in design, specifications and models. The only warranty Baxter makes is the express written warranty extended on the sale or rental of its products. Baxter Medical Systems GmbH + Co. KG

Baxter and Hillrom are trademarks of Baxter International, Inc. or its subsidiaries. EMA-CS84-220004 v2 01/2023