

Baxter

TS7000 U14

OR TABLE



TS7000 U14 UNIVERSAL OPERATING TABLE FOR GREATER PATIENT ACCESS

MAKE YOUR OPERATING ROOM MORE FLEXIBLE

- All the advantages of the TS7000 OR Table but with a shortened back section, allowing greater access for the anesthesia provider regardless of the surgery.
- The modular design along with a comprehensive line of accessories and components provides the ability to configure the table for a wide range of applications and offers surgeons greater visibility in applications such as gynecology, urology and orthopedics.
- Simplified configuration and setup with improved coupling point system.
- Radiolucent carbon fiber components are also available for full patient intraoperative imaging.

COMPONENTS & ACCESSORIES FOR VARIOUS APPLICATIONS



Superior access during urology and gynecology procedures with stirrups and URO cut-out.



Simplified patient positioning in urology and gynecology procedures with four-part leg section.



Advanced orthopedic and trauma applications supported by Extension Unit with MIS hip device.



Spine and neuro procedures on universal OR Table with modular structure and broad range of accessories.

ORDERING INFORMATION	TS7000 U14	
Material-N°	TS7000 U14 (MBW) ¹ Pad TS7000 U14 B	Mat.-No. 2065385 Mat.-No. 2065417
Height Adjustment	59 cm – 114 cm / ± 1 cm (23" – 45" / ± 0.4")	
Trendelenburg/Anti-Trendelenburg	± 45°	
Tilt	± 30°	
Longitudinal Slide	46 cm ± 1 cm (18" / ± 0.4")	
Joint Leg Section Up/Down	+ 90° / - 105°	
Joint Back Section Up/Down	± 90°	
Maximum Load	450 kg (1000 lbs.)	

¹Motorized Base with Wireless Interface

For more information please contact your local Baxter representative.

hillrom.com

Markkinoija Suomessa



Mediq Suomi Oy
Vuoritontunkuja 6, 02200 Espoo
puh. 020 112 1510
asiakaspalvelu@mediq.com, www.mediq.fi

The TS7000 U14 Surgical Table is intended to be used by clinicians and medically qualified personnel. This medical device is a regulated healthcare product which, pursuant to such regulation bears a CE mark. Baxter recommends that you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices. Baxter reserves the right to make changes without notice in design, specifications and models. The only warranty Baxter makes is the express written warranty extended on the sale or rental of its products. Baxter Medical Systems GmbH + Co.KG.

Baxter and Hillrom are trademarks of Baxter International, Inc. or its subsidiaries. EMA-CS85-220003 v2 01/2023