

TruSystem 7500 Hybrid

for Siemens Artis zeego



TruSystem 7500 – free floating carbon table top for TAVI, integrated with Siemens Artis zeego



Therapy and interventional imaging for all clinical situations: The Hybrid OR featuring the TruSystem 7500

Precise

High speed adjustment rates up to 185 mm/sec (7.28 inch/sec) maintaining the isocenter: synchronous movement of the C-arm guarantees a maximal amount of patient safety and efficiency

Flexible

Multi-axial, flexible positioning: various working positions of the C-arm and OR-table provide an increased radiographic area featuring a whole-body coverage that is particularly suitable for angiographic diagnosis

Functional

Precise patient positioning due to special carbon hook coupling points in combination with a broad range of radiolucent equipment: virtually unlimited radiographic coverage during transapical and transfemoral TAVI

Workflow

Enhanced user comfort and efficient working sequences: numerous opportunities and configurations using alternative TruSystem surgical tables with Siemens Artis zeego in parking position

Standard configuration for hybrid operating rooms, main application cardio/vascular

Number	Description	PCS
1854088	Stationary column TruSystem 7500 Hybrid Plus	1
1902265	Floor-mounting plate flat Hybrid	1
1854087	Floor mounting column TruSystem 7500 Hybrid Plus	option
1955207	Installation kit Anker Hybrid ¹	option
1955208	Installation kit Sandwich Hybrid ¹	option
1637573	IR remote control TruSystem 5500	1
1557243	Charging unit mobile AC TruSystem 7500	1
1616204	OR table top Carbon FloatLine V	1
1797654	Pad Carbon FloatLine B	1
1774603	Power Shuttle (H)	1
1774606	Battery Power Shuttle	2
1774608	Charging unit Power Shuttle	1
1489491	Shuttle 3.6 360	option
1272324	Clamp pair for joint plate Carbon	1
1266934	Joint plate Carbon 520	1
4544166	Side rail	2
4544536	Side rail long	1
1939050	Arm Shield, with pad	2
1895203	Carbon Fiber Arm Board 75	2
1697635	Universal arm support 450 T	option
1814876	Catheder tray	1
1946701	Lower body table shield	1
1946702	Wide lift-off top shield 27 cm	1
1946703	Wide lift-off top shield 57 cm	1
1424577	Head section X-RAY wide T	1
1312733	Adapter X-RAY wide	1
1430024	Adapter X-RAY 3-joint wide	option
1708607	Pad for head positioning	1
1708584	Leg positioning pad	1

¹ necessary for the floor mounting column TruSystem 7500 Hybrid Plus

Standard configuration for hybrid operating rooms,
main application neuro, trauma and spinal surgery

Number	Description	PCS
1582739	Operating table top SQ14 X-TRA T	1
1774606	Battery Power Shuttle	2
1774608	Charging unit Power Shuttle	1
1774603	Power Shuttle (H)	1
1489491	Shuttle 3.6 360	option
1430023	Adapter X-RAY 3-joint narrow	1
1317569	Skull clamp X-RAY	1
1317786	Swivel adapter X-RAY	1
1345749	Retainer X-RAY	1
1272324	Clamp pair for joint plate Carbon	1
1279867	Joint plate Carbon 370	1
4544166	Side rail	2
1368880	Coupling piece X-RAY low	1
1300662	Head rest X-RAY	1



SQ14 X-TRA with stationary column TruSystem 7500 Hybrid

Standard configuration for universal purpose

Number	Description	PCS
1846373	OR table top U14 V	option
1846379	OR table top U24 V	option
1846384	OR table top U26 V	option
1851927	Pad back section B	option
1851926	Pad seat section B	option
1789974	Head section 530 double joint V	option
1851979	Pad head section 530 B	option
1789977	Upper back section V	option
1851978	Pad upper back section B	option
1790111	Leg section two parts spread joint V	option
1851976	Pad leg section two parts B	option
1790028	Table top segment Carbon 600 V	option
1770133	Pad TTS Carbon 600 H B	option
1790037	Table top segment Carbon narrow V	option
1873465	Pad TTS Carbon narrow B	option
1489491	Shuttle 3.6 360	option
1774603	Power Shuttle (H)	option
1774606	Battery Power Shuttle	option
1774608	Charging unit Power Shuttle	option

Not all products/ options are available in all countries. For further information about our products or services, please contact your local Trumpf Medical representative or visit our webpage:

www.trumpfmedical.com



Mediq Suomi Oy
PL 115, 02201 Espoo
Puhelin 020 112 1510
asiakaspalvelu@mediq.com, www.mediq.fi

Trumpf Medical reserves the right to make changes without notice in design, specifications and models. The only warranty Trumpf Medical makes is the express written warranty extended on the sale or rental of its products. ©2015 TRUMPF Medizin Systeme GmbH + Co. KG, ALL RIGHTS RESERVED.

TRUMPF Medizin Systeme GmbH + Co. KG
Benzstraße 26 · 82178 Puchheim · Germany · Phone +49 89 80907-0 · E-Mail info@trumpfmedical.com

Flyer ID: 2057078, 30/09/2015

Enhancing outcomes for
patients and their caregivers:

